

BTS Mar165 CCD Pump Down Procedure V2

Subject:

This document describes the steps to pump down the Mar165 CCD head vacuum chamber.

Revision History:

Original 9/27/05 (t madden)

V2: 06/25/200/8 updated photos, added comments (t lutes)

Review Period:

Three years or as needed.

Related Documents:

marusa_marccd_manual.pdf (complete manual from MarUSA)

www.mar-usa.com

BTS Mar165 CCD Pump Down Procedure V2

Tools Needed:

- Mar165 Vacuum Port Adapter (Figure 1)
- Mar165 Vacuum Valve Key (Figure 2)
- Vacuum Pump (Figure 3)
- Mar165 detector and computer running *marccd* GUI

The Mar165 should be pumped down if the pressure reading exceeds 1.0 Torr when the detector is at room temperature. The detector will not cool to the desired temperature, ~ -70C, if the pressure is above this level, and it may take longer to cool if the pressure is 1.0 - 0.5 Torr. The vacuum decreases approximately 0.1 Torr per week. It is possible to pump the chamber to .01 Torr, but that will only increase the lifetime of the vacuum by about a week. Usually 0.1 - 0.2 Torr is sufficient.

Directions also exist in the Mar165 manual at the end, under “Procedure for re-evacuating CCD chamber”, but follow instructions for external vacuum port model.

Detailed Pumping Procedure:

- 1.) Turn on Mar165 controller.
- 2.) Run *marccd* control software.
- 3.) Reboot detector (Configure → Detector → Reboot), but do not start cooling.
- 4.) Verify that CCD is at room temperature.
- 5.) Attach pump and vacuum port adapter to the ¼” vacuum port on the rear of the CCD head labeled “VAC PORT”.
- 6.) Turn on the pump.
- 7.) Remove cover on the “VAC VALVE”.
- 8.) Use the key to **slowly** open the valve, while monitoring the pressure reading. Open the valve a few turns so the pressure falls steadily. The pressure should immediately begin to drop as you open the valve. If not, close the valve and check the connections.
- 9.) Wait until pressure has gone down to 0.1 - 0.01 Torr.
- 10.) Close the VAC VALVE, hand-tight.
- 11.) Turn off pump and remove adaptor.
- 12.) Check that pressure is stable.
- 13.) Turn off the Mar165 detector or start the cooling from *marccd* GUI.

BTS Mar165 CCD Pump Down Procedure V2

Figure 1: Swagelok adapter attached to standard Vacuum hose for the Mar165 vacuum port.



Figure 2: Vacuum valve key. This is in the Mar165 toolkit.



Figure 3: Vacuum Pump. Any style pump that reaches 0.01 Torr is sufficient.

