

Reduction of 9IDC (or 15IDD) USAXS/pinSAXS/WAXS data with Nika

Nika version: 1.70

Introduction

This is step by step procedure for data reduction of pinSAXS data from 9IDC combined USAXS/SAXS instrument. This applies to data collected in 2011-03 APS cycle and later with some changes making this code incompatible to early data. This DOES NOT apply to 2011-02 cycle, when data were saved ONLY in Tiff file format. Some of the tools can be used, but it is not so straight forward. You may want to contact me for help with reduction of data from 2011-02 cycle.


Necessary pre-requisite

You should have USAXS data and reduce them first. If you choose to do so, this customized configuration of Nika may look up the parameters from the USAXS data. Further, you need to have installed Irena package so you can merge the data together.

Procedure

1. Nika configuration

Select from Nika "SAS 2D" > "Instrument configuration" > "9IDC or 15IDD USAXS-SAXS-WAXS"



The image shows two overlapping windows. The background window is titled "Instructions_15IDD" and contains the following text:

Instructions for use of 9IDC (15IDD) special configurations

Decide which data you need to reduce. Instructions are setup specific:

pinSAXS

You may be helped by first reducing your USAXS data and process them to SMR waves, but it is not necessary.

1. Select "pinSAXS" checkbox.
2. Push "Set default methods" button to locate data folder and configure Nika settings common to all pinSAXS experiments.
3. Select one of the images from your pinSAXS measurements in the main 2D panel and double click it (or use button "Ave & Display sel. file(s)") to load
3. Use "Set Experiment Settings" button to read & set values from the wavenote of this file (These are the most likely values for your experiment)
4. Optional: Verify the parameters (Beam center & calibration) using Ag Behenate measurements collected with your data & your notes (some NX file info could be stale)
5. Create mask to mask off the bad points and top few lines on detector using button: "Create bad pix mask"
6. Optional : Either read slit length (button:"Set Slit Length") from one of USAXS measurements and select "Create Smeared Data" (resulting data will be "_usx") - OR - unselect the "Create Smeared Data" (resulting data will be "_270_30"). Likely choose "Delete temp Data" if you are creating slit smeared data.
7. Find "Empty" (aka: Blank) file you want to use and load it in Nika ("Emp/Dk" tab). To process pinSAXS select the 2D data in the main panel and push button "Convert sel. files 1 at time"
8. Select data sets and reduce using "Convert sel. files 1 at time". This should create reasonable data for pinSAXS if you have small-angle scattering. If you have diffraction peaks, select in "Sectors" tab checkbox "Max num points" and if you need it, select "d?" checkbox etc.

Merge "_usx" data using "Data manipulation I" tool from Irena package with USAXS slit smeared data or "_270_30" data with USAXS dsmeared data.

WAXS

1. Select "WAXS" checkbox.
2. Push "Set default methods" button to locate data folder and configure Nika settings common to all WAXS experiments.
3. Select one of the images from your WAXS measurements in the main 2D panel and double click it (or use button "Ave & Display sel. file(s)") to load
3. Use "Set Experiment Settings" button to read & set values from the wavenote of this file (These are the most likely values for your experiment)
4. Select data sets and reduce using "Convert sel. files 1 at time". This should create reasonable data for WAXS. There is no need to empty or other corrections in this case. You may need to modify the "Sectors" tab settings depending if you want the output in Q or in d (or in Theta).


The foreground window is titled "NI1_15IDDConfigPanel" and shows the configuration interface for the 9ID-C (or 15IDD) Nexus file configuration. It includes checkboxes for pinSAXS, WAXS, and 15ID SAXS, buttons for "Set default methods", "Set Experiment Settings", "Set Slit Length", and "Create default pinSAXS mask", and options for "Create Smeared Data" and "Delete temp Data".

Note the instructions next to the panel. Read them if you are nto sure what to do.

You should be able to open this manual when you push button "Open 15IDD manual". If this is not available on your computer, contact me.

Basic Configuration





Use button "Set default methods" which will set part of the parameters, which is independent on the beamline setup. This will force dialog in which you need to find your area detector data. NOTE: the folder with your area detector data was named same as your spec file with `_saxs` at the end.




This will set the basic methods to be used...

Version 2.0 from Sunday, October 25, 2015


See following images. Note, that there may be some changes done to defaults and so trust the code how it sets up the Nika, these images may be obsolete if we decide to make changes to data reduction procedures...






Now select ANY image file and use button “Ave & Display sel. File(s)” on main panel to load it in. Really does not matter which, just from the right folder with the data. The assumption here is, that data collected while using one spec file will have same metadata.


Note two windows, one is image and the other is metadata which were included in that image:



The image may be more informative, if you select “Log Int. display” at the bottom of main Nika panel.



Now use button “Set Experiment Settings” from the “9ID-C (or 15IDD) Nexus File Configuration” to load in default values stored in the wave note.



Some of these are guaranteed to be correct (wavelength, pixel size) and some are likely just close (beam center, distance). We will refine those next.


Beam center/distance fine tuning

Stored beam center and distance between sample and area detector are probably close enough – experience is that while not perfect, any variability has minimum impact on the data.

Therefore I suggest using simply the values stored in the nexus file. If you need more precision or are interested in, there is summary how to refine center and distance at the end of this document.


Further configuration

Set the slit length by selecting “Set Slit Length” button and selecting any USAXS data from the spec file with the same name as your 2D data folder (slit length is same for all data in one USAXS geometry, therefore nearly reliably for all samples in one USAXS spec file).



Slit length should be now set to sensible number (around 0.02 – 0.04 most likely).

You should push the button “Create default pinSAXS mask” for pinSAXS and check if it is appropriate. Most likely is, if not, you need to create your own mask. Check Nika manual how to do that.



This button will create mask covering the horizon tungsten bar we use instead of beamstop and few dead pixels on our detector. It will also select the mask to be used in main panel.



Note, this will overwrite any existing mask in the Igor experiment and so if you want to use other mask (if necessary), then you need to create your own mask and mask yourself points [86,17] and [58,112] as part of that mask.

Select what data to use...

Unless you are real expert, use default values here. Nika will create fake “slit smeared” data which can be merged directly with SMR data from Indra (USAXS). Description is in the back of this file...

Data reduction part:

Generally, if you know what is your empty (aka: Blank) file and have no issues, you can ignore the “Force ...” checkboxes, select the right Empty file in Em/Dk tab and go....




contain the reduced data.

If you have many different Empty/Blanks and need to keep track that you are using the same in pinSAXS as USAXS, you can check “Force ..” checkboxes and then text located at the end of these instructions applies...

Once you have the right Empty/Blank all is left is to select one or more files and click button **“Convert sel. Files 1 at time”**, which will process all selected files, each individually. It will create new folder pinSAXS in the current Igor experiment (if it does not exist) and in this folder will be folder named by the image loaded. Inside will be q_, r_, s_, and w_ waves (“qrs” system), which

Example of data reduced using this procedure...




Merging the data together


From version 2.58 of Irena you should use “Merge data” tool in irena, which follows the same logic, but is simpler.

Note, that if you will follow this routine, some features of the Data Manipulation I tool have been improved by adding “merge” and “Merge 2” buttons, which make this process easier.


To merge the data together we can use the `_usx` data from PinSAXS and SMR data from USAXS, we can use the “Data manipulation I” from Irena.




Select “Indra 2 data” for First data set and “QRS” for second data set, select the same sample measured and “Add data and Graph” button.



Zoom to range where data overlap and subtract background from data set 1 (Modify data 1 area) to match the slope of the curves in the overlap region:




Do not over subtract or the slopes will not match. Select with cursors are where the data overlap well and use "AutoScale" button to scale the data together:




Now you should remove the pinSAXS data for first few points (3-5 points typically) where there are edge effects and other artifacts in pinSAXS data. Use rounded cursor and "Rem Q<CSR(A)" button. Remove

USAXS data at high Q where the noise or background are too large using rectangular cursor and “Rem Q>Csr(B)” button. Make sure the cursors are on the right data set.


Then select “Combine data” as processing method and hit button “Go”.



When happy with the result modify the name of the “New data folder” (may be add behind the name: “_comb”) and save the data.




New data were created. These data have absolute calibration from USAXS instrument. They also are still slit-smearred with the original slit length and can be desmeared or modeled by Irena as slit smearred:



WAXS data processing

Follow same procedure as above, but only two buttons:



Other controls do not apply.

In this case you only need to :

1. Push “Set default methods” and locate the data folder (it will be named by your spec file name with `_waxs` attached).
2. Load any data set
3. Push button “Set Experiment settings”

If you need to subtract background, talk to ilavsky@aps.anl.gov as that requires bit more complicated procedure.


Beam center/distance fine tuning

Stored beam center and distance between sample and area detector are probably close enough – experience is that while not perfect, any variability has minimum impact on the data.


Therefore I suggest using simply the values stored in the nexus file. If you need more precision or are interested in, there is summary how to refine center and distance at the end of this document.

We can refine them next using measurement of Ag Behenate. Note, that you may have more than one Ag Behenate measurement if you had multiple “batches” of samples collected. These are likely close together, so you may choose to use only one setting for all samples, but correct procedure is to optimize Beam center/distance for every batch independently.

Select “Beam center and Geometry cor.” From “SAS 2D” menu.




Find appropriate Ag Behenate data set and hit “Make Image”, select “Log Image?” checkbox in the middle left of the panel:




Note the red dot and circle at the top is position of the beam center stored in the metadata by beamline scientist.

Select tab “Calibrant” and you should see:




This shows that Ag Behenate is predefined as calibrant. The grey lines of Ag Behenate should be between the two red lines around the white line. The white line is estimate of position of the diffraction lines using the distance and beam center from the metadata, two red lines indicate how wide will Nika search for the peak positions. This should be more or less correct. If higher orders are too weak, you may choose to uncheck higher orders of the diffractions here (d5, d4 etc.).

Select tab "Refinement" and select "Refine beam center" and "Refine Sa-Det distance". Run first the refinement ("Run refinement" button) with "Gauss" selected as "Peak shape function". If the fitting to peak profiles starts failing, select "Gausswithslopedback". Gauss should be more stable, but when background is significant, the "Gausswithslopedback" may be necessary.



After fitting, history are will contain message with achieved chi-squared. Now select checkbox “Refine tilts” and run again. The chi-squared should improve somehow, there are likely tilts (about 0.5 deg horizontal, 4.5 deg vertical in my case):



This is done... You now have the best geometry parameters you can have for the batch of data the AG Behenate file is associated with. Close the panel, Sample_Information notebook, and image.


Select what data to use...

Unless you are real expert, use default values here. Nika will create fake “slit smeared” data which can be merged directly with SMR data from Indra (USAXS).

OK, here comes the fun.


1. Area detector data can be reduced as pinhole collimated by doing sector average. Predefined sector is +/-10 degrees around 270 degrees direction (that is vertically down from beam center in Nika notation). These data use as much high-q data as possible reducing the noise, but they are noisy at low-q as number of bins used there is small. They are by default log-binned in Q, further reducing the number of points and therefore noise. NOTE: if you have diffraction lines and want to get high-q resolution, you need to select “Max num points” on the “Sectors” tab to avoid log-binning in q.
2. Area detector data can also be reduced in “Slit smeared” configuration by using lineProf tool around the 270 degrees direction over width of about the slit length. This uses lot more points at low-q but much less points at high q. This tool also cannot be used to bin q scale logarithmically.

See the difference in the images when you add the lines using “Display sects/Lines?” on the main panel.




The slit smeared data can be merged with USAXS SMR data directly, the pinhole data can be merged with the desmeared data only...

So here comes the trick. If you choose “Create Smeared data” checkbox on the 15ID-D Nexus file configuration panel, Nika will generate both of the data listed above and merge them together by using the smeared data at q values smaller than $3 \times \text{slit length}$ and pinhole data at q values higher than $3 \times \text{slit length}$. These data should have the best of both Worlds here: smeared data at low q values to enable merging with SMR USAXS data and lowest noise at high q ...



If you select “Delete temp Data” checkbox, the sector and line profile data will be deleted, reducing amount of data created and therefore clutter. On the other hand, if you want to see differences, you can uncheck it and look at the different data created.

“Force Transmission verification” – since not all samples are homogeneous and USAXS beam size is typically larger than pinSAXS one, the transmission for USAXS and pinSAXS MAY be different. However, we measure ONLY USAXS transmission... To enable user to modify (aka: try, fudge,...) modification of the transmission, here is checkbox which will force user dialog for transmission for each sample when the sample is being reduced. The following window will then appear for each sample:



Use when needed. Carefully...


Force checkboxes...

If you have many different Empty/Blanks and need to keep track that you are using the same in pinSAXS as USAXS, you can check “Force ..” checkboxes and then text located at the end of these instructions applies...

Here we take advantage of the fact, that Nika can look up various parameters from USAXS data reduction and therefore can be helpful.

Select file you want to reduce. Then hit “Convert sel. Files 1 at time” button on main panel.


If Nika finds wrong empty (or no empty at all) it will put out dialog:



Nika is here using name of the USAXS blank name and compares with the empty 2D file loaded. It cleans up extension and order number, also can remove _XYZs which is used by some users to indicate different exposure times. Nika will also check, how many similarly named images are there in the current folder and put out the dialog, if the selection is not unique.

Your choices are to load proper image, keep current or abort in case you are really confused by now.

Choose “Yes” to load proper empty file and dialog is presented:



Note, that the most likely 2D data sets, matching the USAXS blank name are presented at the top. Just to make sure, all 2D data sets are presented below.

Select the right 2D data for empty and it is loaded.

Note: if you do not get dialog for empty data set, it means that Nika code has decided that the current empty name matches sufficiently the name of the USAXS blank. There is no guarantee that this is correct, but I think it would become very annoying to add more dialogs here.

Note that in the history area are records of what was found and used:

Loaded following wave: root:Packages:NexusImport:'11aPPP10_069.hdf5':data

Loaded file 11aPPP10_069.hdf5

For sample : 11aPPP10_069.hdf5 has been found USAXS transmission = 0.93608 in folder : root:USAXS:'11_09_gianf':S85_11aPPP10:

This enables to process large number of data sets quickly and efficiently, as the code is keeping empty file selection under control and is getting transmissions automatically.

If the USAXS data cannot be found, the code will present dialog to user asking for transmission values.

Here are some data which users measured and which were processed by using this code: